

List of CSAN® Pronto® Consumables

Quick Reference Guide



To order : Contact HLS Customer Service at 1-866-669-2313
Email : hls-csr@hlstherapeutics.com



Test Strips

HLS # 80301202

Disposable CSAN® Pronto® test strips.

Available format - Box of 50 test strips



Auto-Check Strip

HLS # 80301207

Reusable CSAN® Pronto® Auto-check strip (quality control) to verify the suitability of the CSAN® Pronto® device.

Available format – Single unit only

To order : Contact CSAN® at 1-800-267-2726

Email : csan@hlstherapeutics.com



Optical Swab

HLS # 80301205

Disposable CSAN® Pronto® optical swab (to clean the optical parts of the CSAN® Pronto® device).

Available format – Single unit only



AC Adaptor

HLS # 80301204

CSAN® Pronto® AC power adapter (to power on the CSAN® Pronto® device and connect it into the power outlet).

Available format – Single unit only



HLS Therapeutics®

HLS Therapeutics Inc.
10 Carlson Court, Suite 701
Etobicoke, Ontario
M9W 6L2
www.hlstherapeutics.com

CLOZARIL, CSAN, and CSAN Patient Care Portal + Design are all registered trademarks of Novartis AG. Pronto is a registered trademark of HLS Therapeutics Inc. All rights reserved. © Copyright 2024 HLS Therapeutics Inc.

January 2024

MLR: CA-CP-00033 ENG

